

Automated Solid Phase Extraction of Dioxins and Furans in Soil and Sediments



Figure 1. Instrumentation used for sample preparation. TurboVap® LV for automated sample evaporation. The RapidTrace® for automated SPE.

Introduction

EPA methods 8260 & 8270 cover the analysis of polychlorinated dibenzo-p-dioxins and polychlorinated dibenzofurans in ground & surface waters. This also now includes 1,4-dioxane and PCBs. This application note will focus on the automated SPE of the sample rather than the detection which usually involves GCMS or LCMS. NB The columns and reagents used for soils and sediment analysis differ to ground water extraction.

Soil and Sediment Sample Preparation

- » A 10–20 g dry weight sample of soil or sediment needs to be mixed with several 20 mL volumes of Hexane : Acetone (90:10) to extract the organic compounds from the matrix. The resulting solvent then needs evaporation using a TurboVap to reduce the volume to approximately 1 mL. The Acetone will be lost by this stage so that the Dioxins and Furans in Hexane can be separated from the interferences to allow low detection limits.
- » An ISOLUTE® NH₂ 500 mg/3 mL SPE column part number 470-0050-B is used to capture the interferences and allow the cleaned up Eluent to be collected. The Eluent can be passed through an anhydrous sodium sulphate cartridge to dry it or a layer of anhydrous sodium sulphate can now be incorporated in new SPE cartridges.

RapidTrace Method

All solvent lines are purged and primed with solvent first. A maximum of 8 solvents can be used to run a wide range of methods. The 4 independent waste solvent lines can be separated if needed. The instrument sample rack has two rows of 10 test tubes for holding the samples and fractions. If multiple fractions are collected from the same sample it will reduce the throughput. The sample preparation time for this dual fraction collection method is under 6 minutes and offers an automated rugged and reproducible solution for busy environmental laboratories.

Step	Source	Destination	Volume (mL)	Flow (mL/Sec)
Condition	Hexane	Organic Waste	5	15
Load	Sample	Fraction 1	1	1
Collect	Hexane	Fraction 1	2	1
Purge - Cannula	Hexane	Cannula Waste	2	30

Common Reagent Table for all Methods

Line No.	Reagent Name	SIP Speed (mL/min)
1	Hexane	30

Waste Name	Abbreviation	SIP Speed (mL/min)
Aqueous Waste	Aq W	Air Push = 2mL
Organic Waste	Org W	Air Push Multiplier = 2
Cannula Waste	Cannula	

Ordering Information

Part Number	Description	Quantity
C50000	RapidTrace®+ Workstation 1 mL and 3 mL (10 columns)	1
C103198	TurboVap® LV 100/120V	1
C103199	TurboVap® LV 220/240V	1
470-0050-B	ISOLUTE® NH2 500 mg/3 mL Sample Volume Columns	50

For the latest application notes visit www.biotage.com

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
 Fax: +86 21 2898 6153
cn_order@biotage.com
cn-1-pointsupport@biotage.com

To locate a distributor,
 please visit our website at
www.biotage.com

Part Number: AN904.V.1.

© 2014 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. Product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

For more information visit www.biotage.com.