

GM9

Fast Plasma Glucose Analyser

for Clinical and Research Applications

APPLICATION AREAS

- Metabolic Studies
- Biochemical Research
- Diabetic Research
- Glucose Clamping



MAIN FEATURES

- Plasma, serum and other aqueous solutions or whole blood via Analog collection systems
- Small sample sizes from 2.5 - 10 µl
- Printed results in under 20 seconds
- No sample turbidity or opacity errors
- Simple YES/NO operation
- Fully sterilizable fluid pathways
- Data output facility
- Compact size
- Fully portable version available

PRINCIPLE OF OPERATION

In the presence of molecular oxygen, β-D-glucose is oxidised by the enzyme glucose oxidase (GOD) to gluconic acid and hydrogen peroxide,



Under the conditions of the assay, the rate of oxygen consumption is directly proportional to glucose concentration.

ANALYTICAL PERFORMANCE

	Accuracy	Linearity	Precision (Within Run)
Glucose >	i) Method comparison vs Hexokinase: $y(\text{Analox}) = 0.985x - 0.14 \text{ mmol/L}$, $r = 0.999$, $n = 156$ ii) Method comparison vs Beckman: $y(\text{Analox}) = 1.005x - 0.07 \text{ mmol/L}$, $r = 0.999$, $n = 123$ iii) Method comparison vs YSI: $y(\text{Analox}) = 1.008x - 0.01 \text{ mmol/L}$, $r = 0.999$, $n = 97$	30.0 mmol/L (540 mg/dl) for 10 µl samples; 50.0 mmol/L (900 mg/dl) for 5 µl samples	C.V. of 1.0 % @ 5 mmol/L (plasma) C.V. of 1.4 % @ 10 mmol/L (plasma) C.V. of 0.85 % @ 12 mmol/L (whole blood)

INSTRUMENT SPECIFICATIONS

Method	> Enzymatic oxygen-rate	Statistical Programmes	> Sequential, giving mean, S.D and C.V.
Sensor	> Clark-type amperometric oxygen electrode	Interface	> Serial data port, optional Windows software available
Sensitivity	> 0.1, or 0.01, selectable		
Reaction Temperature	> 30°C	Power	> 100-250V AC, 50-60Hz, 12-15V DC, 60VA
Display	> 32 character backlit LCD	Dimensions	> Width 23cm, (9 ins) x Depth 29cm, (11½ ins) x Height 15cm, 6¼ ins
Printer	> 16 column dot matrix, 1 line/sec	Weight	> 3.8 kg, 8 lb 6 oz Portable Model 5.9 kg, 13 lb