



ISOLUTE SLE+

ISOLUTE SLE+

ISOLUTE SLE+

ISOLUTE SLE+

PART NO. 820-0140-CG



LOT No. 11108110CC



Content	Size	Content	Reference
ISOLUTE SLE+	50	50	50

ISOLUTE®



VacMaster

VacMaster - 96

# Manual Sample Processing Products and Accessories

---

# Biotage® PRESSURE+ 48 and PRESSURE+ 96

## Positive Pressure Manifolds

### Biotage® PRESSURE+ 48 and PRESSURE+ 96

Biotage® PRESSURE+ manifolds deliver positive pressure, parallel processing for 96 well plates, 1 mL, 3 mL and 6 mL column formats. The systems utilize a consistent, uniform flow of positive pressure to move both low and high viscosity liquids through SPE plates and columns.



### Biotage® PRESSURE+ 96

Part Number	Description
PPM-96	PRESSURE+ 96 Positive Pressure Manifold 96 position. Includes 1 x 1 mL 96 well collection plate, 1 x 2 mL 96 well collection plate, 1 x 10 mL 24 well collection plate (waste plate), 96-column sealing gasket, gas supply adaptor kit (6' of 1/8" i.d. polyethylene tubing and 1/8" and 1/4" NPT connectors), User Manual CD-ROM.

#### Accessories

PPM-A96-CH	PRESSURE+ 96 Tabless 1 mL Column Holder
PPM-A96-SPCR	PRESSURE+ 96 Spacer (16 mm) for $\mu$ Elution Plates
PPM-A96-GSKT	PRESSURE+ 96 Sealing Gasket 96 position
PPM-A96-1024	PRESSURE+ 96 Collection Tray 10 mL 24 well
PPM-GA	PRESSURE+ Gas Supply Adaptor for all models

### Biotage® PRESSURE+ 48

Part Number	Description
PPM-48	PRESSURE+ 48 Positive Pressure Manifold 48 position. Includes Sealing gasket 48 position, gas supply adaptor kit (6' of 1/8" i.d. polyethylene tubing and 1/8" and 1/4" NPT connectors), waste bin rack, waste bin inserts, User Manual CD-ROM. Racks must be ordered separately.

#### Racks

PPM-A48-1RCK	PRESSURE+ 48 SPE Column Rack 1 mL
PPM-A48-3RCK	PRESSURE+ 48 SPE Column Rack 3 mL
PPM-A48-6RCK	PRESSURE+ 48 SPE Column Rack 6 mL
PPM-A48-1232	PRESSURE+ 48 Sample Vial Rack 12 x 32 mm
PPM-A48-1275	PRESSURE+ 48 Collection Rack 12 X 75 mm
PPM-A48-13100	PRESSURE+ 48 Collection Rack 13 x 100 mm
PPM-A48-16100	PRESSURE+ 48 Collection Rack 16 x 100 mm

#### Accessories

PPM-A48-GSKT	PRESSURE+ 48 Sealing Gasket 48 position
PPM-A48-WST	PRESSURE+ 48 Waste Bin Inserts, pack of 3
PPM-A48-WSTRCK	PRESSURE+ 48 Waste Bin Rack with Inserts, pack of 3
PPM-GA	PRESSURE+ Gas Supply Adaptor for all models

### Biotage® VacMaster™ 10 and 20

10 or 20 position vacuum manifolds for SPE, SLE and filtration applications.



#### VacMaster 10 Sample Processing Manifold

Part Number	Description
121-1010	VacMaster-10 Sample Processing Manifold (with 10 mm rack)
121-1012	VacMaster-10 Sample Processing Manifold (with 12 mm rack)
121-1016	VacMaster-10 Sample Processing Manifold (with 16 mm rack)
121-1025	VacMaster-10 Sample Processing Station (with 25 mm rack)
121-1027	VacMaster-10 Sample Processing Station (with 27 mm rack)

#### VacMaster 10 Replacement Parts

Part Number	Description	Qty.
121-1039	Replacement tank with fittings	1
121-1045	Replacement lid with fittings	1
121-1030	Silicone Lid Gasket	1
121-1133	Replacement rack, 12 mm	1
121-1134	Replacement rack, 16 mm	1

#### VacMaster 20 Sample Processing Manifold

Part Number	Description
121-2010	Sample Processing Manifold (with 10 mm rack)
121-2012	Sample Processing Manifold (with 12 mm rack)
121-2016	Sample Processing Manifold (with 16 mm rack)
121-2020	VacMaster-20 Sample Processing Station (with 20 mm rack)

#### VacMaster 20 Replacement Parts


Part Number	Description	Qty.
121-2068	Replacement tank with fittings	1
121-2075	Replacement lid with fittings	1
121-2059	Silicone Lid Gasket	1
121-2161	Replacement rack, 10 mm	1
121-2162	Replacement rack, 12 mm	1
121-2163	Replacement rack, 16 mm	1

#### VacMaster Accessories

Part Number	Description	Qty.
121-0010	PTFE 'T' valve (for connecting 2 tanks)	1

### PTFE Stopcock Option and Spare Parts

Part Number	Description	Qty.
1 121-0009	Universal PTFE stopcock	10
2 121-0001	PTFE stopcock/needle unit	10
3 121-0002	PTFE needle	10
4 121-0003	Stainless steel needle	20
121-0004	Stainless steel needle retainer	10
121-0005	Port sealing plugs	30
121-0009-S	PTFE stopcock - pressure positive	10



# Biotage® VacMaster™ 96

## Sample Processing Manifold

### Biotage® VacMaster™ 96

Vacuum manifold for processing 96-well SPE, SLE or filtration plates.

Part Number	Description	Qty.
121-9600	VacMaster-96 Sample Processing Manifold (without vacuum control)	1

#### Vacuum Control Options

Part Number	Description	Qty.
121-9601	VacMaster VCU-1 Vacuum Control Unit	1
121-9602	VacMaster VCU-2 Vacuum Control and Generation Unit	1

#### 96 Well Collection Plates-Deep Well

Part Number	Description	Qty.
121-5202	Collection plate, 1 mL, square	50
121-5203	Collection plate, 2 mL, square	50
121-5213	Collection Plate, 2 mL, round	50



### Replacement Parts


Part Number	Description	Qty.
121-9612	VacMaster-96 replacement gasket	1
121-9613	VacMaster-96 replacement o-ring	1
121-9614	Collection plate spacer (2 mm) for deep well collection plates	1
121-9610	Array insert (6 mm), acetal	1
121-9615	Collection plate spacer (29 mm) for shallow well collection plates	1
121-9611	VacMaster-96 insert (12 mm) for "shallow skirt" plates	1



## Biotage® VacMaster™ Drying Adaptor

Connect to laboratory air or nitrogen supply to dry 10 or 20 SPE columns simultaneously.

Part Number	Description	Qty.
124-1001	VacMaster-10 Drying Adaptor for 1, 3 & 6 mL columns	1
124-2001	VacMaster-20 Drying Adaptor for 1, 3 & 6 mL columns	1
124-2002	VacMaster-20 Drying Adaptor for 15 and 20 mL columns	1



## Biotage® VacMaster™ Trap Kit

Waste traps should be installed between the outlet of the VacMaster™ sample processing manifold and the vacuum source, trapping all waste liquids exiting the manifold. Compatible with VacMaster 10 & 20 and -96 processing manifolds, VacMaster Trap Kits are available with 1 L or 10 L capacity.

Part Number	Description	Qty.
121-2095	VacMaster Trap Kit 1 L	1
121-2195	VacMaster Trap Kit 10 L	1



## Biotage® VacMaster™ Large Volume Extraction (LVE) Kit

For unattended loading of large volume samples. Inert PTFE tubing prevents sample contamination.

Part Number	Description	Qty.
121-2090	VacMaster LVE kit for 1, 3 and 6 mL columns	1
121-2091	VacMaster LVE Kit for 15 mL and 25 mL SPE columns	1
121-2092	VacMaster LVE Kit for 70 mL SPE columns	1
121-2094	VacMaster LVE Kit XL	1



# Biotage® Gravity Rack

## Biotage® Gravity Rack

Process up to 20 samples simultaneously under gravity.

Part Number	Description	Qty.
123-2016	Gravity Rack with 16 mm collection tube rack	1
123-2019	Gravity Rack with 19 mm collection tube rack	1

### Accessories

Part Number	Description	Qty.
121-0009	Universal PTFE stopcock	10
121-0001	PTFE stopcock/needle unit	10
121-0002	PTFE needle unit	10
121-0003	Stainless steel needle	20
121-0004	Stainless steel needle retainer	10

