

ISOLUTE® QuEChERS

Fast and Efficient Clean up of
Complex Food Samples

ISOLUTE® QuEChERS

ISOLUTE® QuEChERS

ISOLUTE® QuEChERS (**Q**uick, **E**asy, **C**heap, **E**ffective, **R**ugged and **S**afe) products provide simple clean up of complex samples using salt assisted extraction and partitioning, followed by dispersive SPE. The pre-weighed extraction and clean up tubes conform to AOAC and EN methodologies, and include options for waxed and highly pigmented fruit and vegetable samples.



Extraction Tubes

Part Number	Description	Tube	Pack Size	MgSO ₄	Na Acetate	Na Citrate	Na Citrate sesquihydrate	NaCl
Q0010-15V	15 g QuEChERS AOAC 15 mL Extraction Tube	15 mL	25	6 g	1.5 g			
Q0020-15V	10 g QuEChERS EN 15 mL Extraction Tube	15 mL	25	4 g		1 g	0.5 g	1 g

Clean up Tubes

Part Number	Description	Tube	Pack Size	PSA	MgSO ₄ (purest)	C18(EC)	GCB
Q0030-15V	AOAC Fruit and Vegetables Clean up Tube	15 mL	25	400 mg	1200 mg		
Q0035-15V	EN Fruit and Vegetable Clean up Tube	15 mL	25	150 mg	900 mg		
Q0050-15V	AOAC Waxed Fruit and Vegetables Clean up Tube	15 mL	25	400 mg	1200 mg	400 mg	
Q0060-15V	EN Waxed Fruit and Vegetables Clean up Tube	15 mL	25	150 mg	900 mg	150 mg	
Q0070-15V	AOAC Pigmented Fruit and Vegetables Clean up Tube	15 mL	25	400 mg	1200 mg		400 mg
Q0080-15V	EN Pigmented Fruit and Vegetables Clean up Tube	15 mL	25	150 mg	900 mg		15 mg
Q0090-15V	EN Highly Pigmented Fruit and Vegetables Clean up Tube	15 mL	25	150 mg	900 mg		45 mg

Centrifuge Tubes

Part Number	Description	Tube	Pack Size
Q0000-50V	50 mL Centrifuge Tube with Rack	50 mL	25