



PermeGear is pleased to announce our new 100 μ l to 1000 μ l Pipettes (Part # PP-01) & Pipette Tip Assemblies. They are a three piece flexible assembly, the first designed specifically for insertion into the center of Franz Cells for more accurate sampling. The tip assemblies are available in packages of 100 (Part # PT-100) or 250 (Part# PT-250).

Because of this unique - never before introduced, slender design, you will need to dial in the pipette at a slightly higher volume to achieve your target aliquot. Once a PermeGear Pipette is dialed in, it will maintain $\pm 2\%$ accuracy.

Please keep in mind many factors can affect pipetting results, i.e. temperature, humidity, angle, rate at which sample is pulled, operator techniques etc... Remember to sample slowly for micropipetting.

PermeGear Pipettes:

- Light weight
- Soft spring system for extremely low pipetting forces
- Controlled volume setting to prevent accidental volume changes
- Large volume display
- Contoured shape: fits either large or small hands
- Large pushbutton, rounded and freely rotating
- Finger hook: takes the weight, for a more relaxed grip
- 3-position tip ejector button for left & right handed users
- Easy on site calibration (calibration key included)
- A unique serial number
- UV resistance even prolonged exposure
- Fully autoclavable (120°C/0.1MPa/20 min)
- Calibration key included in the box



If you need additional information or have any questions, please e-mail us at: support@permegear.com