

# Biotage<sup>®</sup> PRESSURE+ 48 and PRESSURE+ 96

---

Positive Pressure Manifolds

# Biotage® PRESSURE+ 48 and PRESSURE+ 96

## Positive Pressure Manifolds

### Biotage® PRESSURE+ 48 and PRESSURE+ 96

Biotage® PRESSURE+ manifolds deliver positive pressure, parallel processing for 96 well plates, 1 mL, 3 mL and 6 mL column formats. The systems utilize a consistent, uniform flow of positive pressure to move both low and high viscosity liquids through SPE plates and columns.



### Biotage® PRESSURE+ 96

Part Number	Description
PPM-96	PRESSURE+ 96 Positive Pressure Manifold 96 position. Includes 1 x 1 mL 96 well collection plate, 1 x 2 mL 96 well collection plate, 1 x 10 mL 24 well collection plate (waste plate), 96-column sealing gasket, gas supply adaptor kit (6' of 1/8" i.d. polyethylene tubing and 1/8" and 1/4" NPT connectors), User Manual CD-ROM.

#### Accessories

PPM-A96-CH	PRESSURE+ 96 Tabless 1 mL Column Holder
PPM-A96-SPCR	PRESSURE+ 96 Spacer (16 mm) for $\mu$ Elution Plates
PPM-A96-GSKT	PRESSURE+ 96 Sealing Gasket 96 Position
PPM-A96-1024	PRESSURE+ 96 Collection Tray 10 mL 24 Well
PPM-GA	PRESSURE+ Gas Supply Adaptor for all Models

### Biotage® PRESSURE+ 48

Part Number	Description
PPM-48	PRESSURE+ 48 Positive Pressure Manifold 48 Position. Includes sealing gasket 48 position, gas supply adaptor kit (6' of 1/8" i.d. polyethylene tubing and 1/8" and 1/4" NPT connectors), waste bin rack, waste bin inserts, User Manual CD-ROM. Racks must be ordered separately.

#### Racks

PPM-A48-1RCK	PRESSURE+ 48 SPE Column Rack 1 mL
PPM-A48-3RCK	PRESSURE+ 48 SPE Column Rack 3 mL
PPM-A48-6RCK	PRESSURE+ 48 SPE Column Rack 6 mL
PPM-A48-1232	PRESSURE+ 48 Sample Vial Rack 12 x 32 mm
PPM-A48-1275	PRESSURE+ 48 Collection Rack 12 x 75 mm
PPM-A48-13100	PRESSURE+ 48 Collection Rack 13 x 100 mm
PPM-A48-16100	PRESSURE+ 48 Collection Rack 16 x 100 mm
415690SP	TurboVap P+ Multi Rack (48 Positions, 12–13 mm Tubes)
415675SP	TurboVap P+ Rack (48 Positions, 16 x 100 mm Tubes)

#### Accessories

PPM-A48-GSKT	PRESSURE+ 48 Sealing Gasket 48 position
PPM-A48-WST	PRESSURE+ 48 Waste Bin Inserts, pack of 3
PPM-A48-WSTRCK	PRESSURE+ 48 Waste Bin Rack with Inserts, pack of 3
PPM-GA	PRESSURE+ Gas Supply Adaptor for all Models